











Length	Type	
32mm	<b>10uL standard</b>	
46mm	<b>10uL extra long (bevelled)</b>	
46mm	<b>10uL extra long (thin &amp; sharp)</b>	
51mm	<b>200uL yellow, bevelled</b>	 standard 200uL, normally used 200uL in bulk pack, and 20uL with filter as filter tips
51mm	<b>200uL clear, bevelled</b>	 bigger body, gives 200uL volume even with filter, and 50uL/100uL filter tips are made of this type too.
60mm	<b>200uL extra long (straight)</b>	
90mm	<b>200uL super long</b>	 For some specific application
78mm	<b>1000uL standard blue</b>	
88mm	<b>1000uL (medium)</b>	 For some specific application
102mm	<b>1,000uL extra long (so called, 1,250uL) blue</b>	 Long enough to give 1,250uL max volume

**10uL, DNase/RNase free, 32mm**



Part NO.	Volume	With Filter	Sterile	Package Type	Pack
UFPT-10R	0.1-10 ul	-	Yes	rack pack	96/rack, 50racks/case
UFPT-F-10	0.1-10 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-REF-10	0.1-10 ul	-	Yes	refill pack	96/rack, 10racks/unit, 10units/case
UFPT-10B	0.1-10 ul	-	-	bulk pack	1,000/bag, 40bags/case
UFPT-F-10B	0.1-10 ul	Yes	-	bulk pack	1,000/bag, 40bags/case

**10uL extra long design, DNase/RNase free, 46mm**



Part NO.	Volume	With Filter	Sterile	Package Type	Pack
UFPT-10ER	0.1-20 ul	-	Yes	rack pack	96/rack, 50racks/case
UFPT-F-10E	0.1-20 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-10EB	0.1-20 ul	-	-	bulk pack	1,000/bag, 40bags/case
UFPT-F-10EB	0.1-20 ul	Yes	-	bulk pack	1,000/bag, 40bags/case

**10uL extra long, thin and sharp design, DNase/RNase free, 46mm**



Part NO.	Volume	With Filter	Sterile	Package Type	Pack
UFPT-10ETR	0.1-20 ul	-	Yes	rack pack	96/rack, 50racks/case
UFPT-F-10ET	0.1-20 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-REF-10ET	0.1-20 ul	-	Yes	refill pack	96/rack, 9racks/unit, 10units/case
UFPT-10ETB	0.1-20 ul	-	-	bulk pack	1,000/bag, 40bags/case
UFPT-F-10ETB	0.1-20 ul	Yes	-	bulk pack	1,000/bag, 40bags/case

**200uL standard, bevelled, DNase/RNase free, yellow, 51mm**



Part NO.	Volume	With Filter	Sterile	Package Type	Pack
UFPT-200R	1-200 ul	-	Yes	rack pack	96/rack, 50racks/case
UFPT-F-20	1-20 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-REF-200	1-200 ul	-	Yes	refill pack	96/rack, 12racks/unit, 10units/case
UFPT-200B	1-200 ul	-	-	bulk pack	1,000/bag, 20bags/case
UFPT-F-20B	1-20 ul	Yes	-	bulk pack	1,000/bag, 20bags/case

**200uL standard wide, Bevelled, DNase/RNase free, 51mm**

Part NO.	Volume	With Filter	Sterile	Package Type	Pack
UFPT-F-50	1-50 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-F-100	1-100 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-F-200	1-200 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-F-50B	1-50 ul	Yes	-	bulk pack	1,000/bag, 20bags/case
UFPT-F-100B	1-100 ul	Yes	-	bulk pack	1,000/bag, 20bags/case
UFPT-F-200B	1-200 ul	Yes	-	bulk pack	1,000/bag, 20bags/case

**200uL extra long or 300uL, straight, DNase/RNase free, 60mm**

Part NO.	Volume	With Filter	Sterile	Package Type	Pack
UFPT-200ER	1-300 ul	-	Yes	rack pack	96/rack, 50racks/case
UFPT-F-100E	1-150 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-F-200E	1-300 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-REF-200E	1-300 ul	-	Yes	refill pack	96/rack, 10racks/unit, 10units/case
UFPT-200EB	1-300 ul	-	-	bulk pack	1,000/bag, 20bags/case
UFPT-F-200EB	1-300 ul	Yes	-	bulk pack	1,000/bag, 20bags/case

**1000uL standard, DNase/RNase free, blue, 78mm**

Part NO.	Volume	With Filter	Sterile	Package Type	Pack
UFPT-1000R	100-1,000 ul	-	Yes	rack pack	96/rack, 50racks/case
UFPT-F-1000	100-1,000 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-1000B	100-1,000 ul	-	-	bulk pack	1,000/bag, 10bags/case
UFPT-F-1000B	100-1,000 ul	Yes	-	bulk pack	1,000/bag, 10bags/case

**1000uL extra long or 1250uL, DNase/RNase free, blue, 102mm**

Part NO.	Volume	With Filter	Sterile	Package Type	Pack
UFPT-1000ER	100-1,200 ul	-	Yes	rack pack	96/rack, 50racks/case
UFPT-F-1000E	100-1,200 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-REF-1000E	100-1,200 ul	-	Yes	refill pack	96/rack, 8racks/unit, 10units/case
UFPT-1000EB	100-1,200 ul	-	-	bulk pack	1,000/bag, 10bags/case
UFPT-F-1000EB	100-1,200 ul	Yes	-	bulk pack	1,000/bag, 10bags/case

**Some specific application****Super long 200uL, DNase/RNase free, 90mm**

Part NO.	Volume	With Filter	Sterile	Package Type	Pack
UFPT-200SR	1-200 ul	-	Yes	rack pack	96/rack, 50racks/case
UFPT-F-200S	1-200 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-200SB	1-200 ul	-	-	bulk pack	1,000/bag, 10bags/case
UFPT-F-200SB	1-200 ul	Yes	-	bulk pack	1,000/bag, 10bags/case

**Medium long 1000uL, DNase/RNase free, 88mm**

Part NO.	Volume	With Filter	Sterile	Package Type	Pack
UFPT-1000MR	100-1,000 ul	-	Yes	rack pack	96/rack, 50racks/case
UFPT-F-1000M	100-1,000 ul	Yes	Yes	rack pack	96/rack, 50racks/case
UFPT-1000MB	100-1,000 ul	-	-	bulk pack	1,000/bag, 10bags/case
UFPT-F-1000MB	100-1,000 ul	Yes	-	bulk pack	1,000/bag, 10bags/case

