

**Cryogenic vials, self-standing, sterile**


Part NO.	Volume	Thread Type	Pack
STSS-512IN	2.0ml	Internal	50/bag, 20bags/case
STSS-512EX	2.0ml	External	50/bag, 20bags/case
STSS-513EX	5.0ml	External	50/bag, 20bags/case

**Cryogenic vials, self-standing, bottom barcode, sterile**

Part NO.	Volume	Thread Type	Pack
STSS-514EX	0.5ml	External	50/bag, 20bags/case
STSS-515EX	1.0ml	External	50/bag, 20bags/case
STSS-516EX	1.5ml	External	50/bag, 20bags/case
STSS-517IN	2.0ml	Internal	50/bag, 20bags/case
STSS-517EX	2.0ml	External	50/bag, 20bags/case
STSS-518EX	5.0ml	External	50/bag, 20bags/case

**PC Cryogenic Boxes, individual pack, stable from -196°C to 121°C**

Part NO.	Height	Wells	Dimensions (mm)	Pack
STSS-BOX5213	2in	81	133*133*53mm	30/case
STSS-BOX5214	3.75in	81	133*133*95mm	30/case
STSS-BOX5210	2in	100	133*133*53mm	30/case

**PP Cryogenic Boxes, stable from -90°C to 121°C**

Part NO.	Height	Wells	Dimensions (mm)	Pack
STSS-BOX5212	2in	81	133*133*53mm	30/case
STSS-BOX5218	2in	100	133*133*53mm	30/case

**Cardboard Cryogenic Boxes, stable from -196°C to 121°C**

Part NO.	Height	Wells	Dimensions (mm)	Pack
STSS-BOX5211	2in	81	133*133*53mm	30/case